**H**

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### INFRARED LAMP

#### Model 3000

Thank you for purchasing a Momert product. We hope you will be satisfied with our product throughout its service life. Please study the entire Operating Manual carefully before you start using the product. Keep the manual in a safe place for future reference. Make sure other people using the product are familiar with these instructions.

#### Technical parameters:

Voltage: 230V~50Hz  
Power: 100W  
Bulb: Suntown - Xwatt R95 100W

#### Legend:

Attention! Risk of injury or damage to health.

Danger hot surface.

#### APPLICATION:

This infrared lamp is only intended for radiating the human body. Infrared radiation, a range between 780-1400 nm, causes warmth to be transported to the human body. Use – beauty care: where the temperature increase stimulates the blood supply to the skin and therewith releases the pollutants and slough, which results in a softer, silky skin.

However, first ask your family doctor whether application in the individual case is advisable from medical standpoint.

**Attention:** This product was not made for people with reduced physical, detection or mental capabilities, or people without experience or knowledge of the unit (children included), only with supervision or informing of a person who knows the use of product, and responsible for their safety.

In case of use by children the supervision is important, to make sure they don't play with the unit.

Cleaning and user maintenance shall not be made by children without supervision.

#### IMPORTANT INFORMATION! PLEASE RETAIN FOR LATER USE.

Please read the information carefully and keep it accessible for other users.

The lamp housing and the infrared lamp become very hot during operation. There is a risk of burning if these are touched. Always unplug the power cord and let the appliance cool down for at least 20 minutes before touching, wrapping, or covering it.

- Remove all the packaging material before the appliance is used.
- The appliance may only be connected to the mains voltage on the type plate.
- Do not leave the appliance unattended when it is switched on.
- Do not use it on people sensitive to heat. The sensitivity to heat can be reduced or increased in the following cases:
  - diabetic patients;
  - people with sleepiness, dementia or difficulties concentrating;
  - people with illness-related skin changes or people with scarred skin areas in the application area;
  - people with allergies;
  - children and older people;
  - after consuming medications or alcohol.
- The appliance is not to be used by:
  - persons suffering from sunburn;
  - persons under medical care for diseases that involve photosensitivity;
  - persons receiving photosensitising medications;
  - a statement that if unexpected side effects, such as itching, occur within 48h of the first session of using an appliance, medical advice should be sought prior to further exposure.
- Advice should be sought from a physician prior to use of the treatment procedures.
- Always limit the period of use and check the reaction of the skin.
- Medications, cosmetics or food can lead to an over-sensitive or allergic reaction of the skin under certain conditions. Radiation

must be ended immediately in this case.

- If the appliance is operated continuously, caution and attention are advisable – never fall asleep during radiation!
  - Prolonged use can cause burns to the skin.
  - Do not touch the infrared lamp or unscrew it when it is connected to the mains.
  - Children do not understand the risks associated with electrical appliances. Make sure the appliance is not used by children without appropriate supervision.
  - The appliance may be used only under supervision.
  - Inflammable objects must be kept at a distance of at least 1,0 m from the infrared lamp. Do not hang the appliance on the wall or from the ceiling.
- DANGER!**
- **Electricity in water is extremely dangerous! Do not use the appliance in a damp environment (such as a bathroom or shower or swimming pool). Make sure there is no liquid inside or outside the device.**
  - **Do not touch the appliance with wet hands when it is connected to the mains; water should not be sprayed on the appliance. Operate it only when it is completely dry. Never immerse into water or any other liquid.**
  - Protect the appliance against blows and impact.
  - Do not remove the plug from the socket by pulling the cable.
  - In case the cable or the casing is damaged, please contact the customer service, as special tools are required for repairing it.
  - Do not use the appliance outdoors. Only use it for its intended purpose
  - In case of misuse, wrong handling or incorrect repair no liability will be taken for possible damages. Claims under guarantee period are also excluded in such cases.
  - The appliance is not suitable for commercial usage.

#### USE:

Remove all packaging material before using the appliance. Before use, check whether the infrared lamp is properly screwed in. Always place the appliance on a stable, level surface and leave the appliance in a radius of at least 20 cm to prevent overheating. Before using make sure the connecting cable is not damaged. It must not come into contact with the hot electric light bulb. **The distance between the infrared lamp and the radiated part of the body should not be less than 60 cm for models 3000 and 3002, or at least 80 cm for model 3001.** Always keep your eyes closed during face treatment to avoid injury (never look directly into infrared light).

**Protective goggles are recommended, especially for children.** The recommended duration of treatment is maximum 12 minutes.

The lamp should be connected to the mains only during operation. After finish using unplug it and let it cool down before storage. We recommend selecting a shorter radiator period to start with. However, the skin of some people may show an over-sensitive (e.g. with extreme redness, blistering, itching, heavy sweating) or allergic reaction to heat radiation, even when the lamp is used correctly. Examine the radiated section of the body regularly. Discontinue radiation immediately if any signs of an over-sensitive or allergic reaction appear and consult physician.

**Do not use the appliance in a room where the temperature is lower than 10°C or higher than 35°C.**

#### CLEANING/CHANGING THE INFRARED LAMP:

Before every cleaning and when exchanging the infrared lamp, the appliance must be switched off, disconnected from the mains, and let it cool down for 20 minutes before touching.

The appliance can be cleaned with a damp cloth. Do not use any cleaning agents containing solvents on the appliance. Ensure that no water penetrates the appliance. Use the appliance again only when it is completely dry. Please note that only infrared lamps of the original type may be used as replacements.

**Bulbs are not included in warranty.**

#### ENVIRONMENTAL PROTECTION:

- Packaging materials and obsolete appliances should be recycled.
- The transport box may be disposed of as sorted waste.
- Polyethylene bags shall be handed over for recycling.

**Appliance recycling at the end of its service life:**

A symbol on the product or its packaging indicates that this product should not go into household waste. It must be taken to the collection point of an electric and equipment recycling facility. By making sure this product is disposed of properly, you will help prevent the negative effects on the environment and human health that would otherwise result from inappropriate disposal of this product. You can learn more about recycling this product from your local authorities, a household waste disposal service or in the shop where you bought this product.

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